

## Purification of DNA via spin column (for DNA in solution)

### GENERAL NOTES:

This protocol is good for restriction fragments, small digested plasmids, and PCR products. It is NOT good for genomic or chromosomal DNA. Use the FastPrep reagents or CTAB-phenol type prep instead for genomic DNA. Digested plasmids work well up to about 5 kb, but the yield quickly falls off past this point. Generally, plasmids don't purify as well as PCR products on these columns.

<<put on gloves + safety glasses. PB contains guanidine which is somewhat toxic >>

1. Figure out what volume of sample you have. If this is less than 100  $\mu\text{l}$ , then make it up to 100  $\mu\text{l}$  with TE. The volume you now have will be called "1 volume". Add 5 volumes of buffer PB to the DNA solution. e.g. if you have a 100  $\mu\text{l}$  restriction digest, add 500  $\mu\text{l}$  PB.
2. Place a silica-based spin column (e.g. 'Econospin') into its 2 ml catch tube (if it isn't already set up that way). Load up to 750  $\mu\text{l}$  of the DNA-PB mixture onto the column. Spin at  $\sim 10,000$  g for 30 sec. Discard the flow-through into culture waste.
3. If you still have more DNA-PB mixture left, repeat the previous step until all of the mixture has been put thru the column. The columns will hold a total of  $\sim 10$   $\mu\text{g}$  DNA, which is a lot!

<<Can take off glasses and gloves now >>

4. Add 750  $\mu\text{l}$  of buffer PE to the column, spin  $\sim 10,000$  g for 30 sec, discard flow-through.
5. Repeat step 4.
6. Spin again for 30 sec to remove all traces of PE from the column. Discard both the flow-through and catch tube, and transfer the spin column onto a clean Kimwipe. Leave the column lid open. Transfer Kimwipe to 60°C oven, and allow to dry for 10 min.
7. Transfer spin column to a sterile 1.5 ml Eppi tube, and add 20-50  $\mu\text{l}^*$  of EB buffer to the centre of the spin column – ie on the membrane, not the walls of tube. Allow to sit for 2 min. Spin at  $\sim 10,000$  g for 1 min, retain Eppi tube with DNA solution in EB, discard spin column.

\* Usually we want high concentration rather than high yield, so use 20  $\mu\text{l}$ . If max. yield is important, or if you have lots of DNA, use 50  $\mu\text{l}$ . Note that you lose approx 5  $\mu\text{l}$  EB during the procedure.

### DNA PURIFICATION SOLUTIONS (SPIN COLUMN) – based on recipes at [openwetware.org](http://openwetware.org)

**Buffer PB:** 5 M guanidine HCl, 30% isopropanol. Not sterilised.

Prepare in a sterile bottle or tube containing the correct amount of autoclaved RO water.

**Buffer PE:** 10 mM Tris-HCl, pH 8, 80% ethanol. Not sterilised.

Add appropriate vol. of autoclaved 1 M Tris-HCl to the correct amount of autoclaved RO water in a sterile bottle or tube, then add 100% ethanol to give 80% final conc.

**Buffer EB:** 10 mM Tris-HCl, pH 8. Autoclaved.

Add appropriate vol. of autoclaved 1 M Tris-HCl to the correct amount of autoclaved RO water in a sterile bottle or tube.