

Preparation of electrocompetent *Mycobacterium smegmatis* & electroporation

Protocol for making electrocompetent cells

1. Inoculate a 5 mL culture of *M. smegmatis* mc²-155 in LB broth containing 0.05% Tween 80 (add after autoclaving from filter-sterile stock) and grow until saturation (~2 days). If there are plasmids in the *M. smeg* cells (e.g. pJV53), add antibiotic too (e.g. Km at 25 µg/ml). You can grow *M. smeg* at either 30°C or 37°C.
2. Subculture into 50 mL of MSM minimal medium in a 250 ml flask (see elsewhere in lab protocols for MSM recipe) containing 1% succinate and 0.05% Tween 80. Start at a final OD₆₀₀ = 0.02 and grow with shaking overnight. If the culture has plasmids in it like pJV53, it will grow more slowly, and a higher initial OD₆₀₀ might be helpful (e.g. 0.1).
3. On the following day, when cells are grown up, but not yet in stationary phase (aim for OD₆₀₀ between 0.6-0.8), transfer the cells to 50 mL Falcon tubes, centrifuge at ~ 4000 g for 10 min at 4°C, discard the supernatant. For cultures containing pJV53, first add acetamide (0.2% w/v) and grow a further 3 hours to induce the recombination enzymes before centrifuging the cells.
4. Add 20 ml ice cold 10% glycerol + 0.05% Tween solution, centrifuge, discard supernatant. Repeat this step two more times. This removes salt, which is bad for electroporation.
5. Resuspend cells in 2 mL glycerol-Tween solution, divide into 200 µl aliquots, freeze at -80°C.

Protocol for electroporation

1. Thaw aliquots of electrocompetent cells (from above) equivalent to half of your number of DNA samples, divide into 100 µl aliquots in Eppi tubes and keep on ice.
2. Add plasmid DNA to the cells. Depending how clean and how concentrated your plasmid is, and what the aim of the experiment is, you may be adding 1-10 µl of plasmid (could be 10-1000 ng, depending on the experiment). If unsure, start with less plasmid and see how you go. Too much plasmid will cause sparking (bad).
3. Allow 10 min on ice for the DNA to adsorb to the cells. While you are waiting, add the electroporation cuvettes to the ice to allow them to chill. We usually use 2 mm gap cuvettes for this. Careful not to get ice inside the cuvettes (ice is not sterile!).
4. Transfer each cell+DNA mix into an electro-cuvette, tap gently to get liquid to bottom. Take samples on ice to the electroporation machine. Set the machine for 2.5 kV, 800 Ω, 25 µF.
5. Wipe down the outside of the cuvette with paper towel, and put it into the holder. Ensure that the metal sides of the cuvette are contacting the electrodes. For some combinations of machine+cuvette, you need to put a spacer in the bottom of the holder to ensure good contact (we use a 5 mm glass bead for this).
6. Hit the button to initiate the electric pulse. Note the time constant – this should be between 15-25 milliseconds. If it sparks, you can try adding 100 µl of 10% glycerol and re-pulsing, but usually this means you need to add less DNA and/or better-quality DNA.
7. Immediately add 1 ml LB-Tween to each cuvette, mix with p1000, transfer to Eppi tube, recover for 2 h with shaking at 30°C or 37°C.
8. Spread-plate 100 µl on the appropriate antibiotic agar plate. Pellet remaining cells, resuspend in a drop of the residual liquid, then spread on a second plate. This gives a mini-dilution series, so hopefully one of the plates will have nice well-separated pickable colonies (this may not matter for some experiments, but does matter for others).